



Appliance Product Information Catalog



When you use our services you can be assured of the quality of our appliances as well as the technical advice and support we provide.

Experience *the* difference...
your patients *will* appreciate *the* results!

Inman Orthodontic Laboratories Inc.
3953 NW 126th Avenue
Coral Springs, FL 33065

www.inmanortho.com
www.inmandligner.com



Welcome to Inman Orthodontic Lab

Our Pledge

Inman Orthodontic Laboratories Inc stands behind every appliance we fabricate.

When you use our services you can be assured of the quality of our appliances as well as the technical advice and support we provide.

I promise to overwhelm you with the personal attention our customer service department will provide you and how easy it is to work with us.

Since 1977 Inman Orthodontic Laboratories Inc has been providing exceptional service to the dental profession and as we have grown we have invested in new technologies. One new technology is our use of Laser Welding, which allows us to fabricate better, stronger appliances by reducing or eliminating the use of silver solder. Our LabNet software allows us instant access to any information you or your staff might need regarding any services we have provided for you.

Our phone system technology virtually eliminates the chance that you will ever experience a busy signal when contacting us, resulting in a quick response to your questions.

As an innovator in the design and construction of Orthodontic Appliances we work hard to keep up with new ideas, trends, and technology.

It is our mission not only to provide unsurpassed quality and customer service but to continue to improve on the design and function of orthodontic appliances.

Sincerely,

Donal P. Inman CDT President
Inman Orthodontic Laboratories Inc.

Experience the difference...
your patients will appreciate the results!



WELCOME

Table of Contents


- 1 Our Pledge
- 2-3 Inman Aligner Family of Appliances
- 4-5 Inman Aligner – *Adjusting the Appliance*
- 6-7 Inman Power Component
- 8-9 Arch Development
- 10 Holding and Pediatric
- 11 Final Retention
- 12 Auxiliary Products – *Clasps and Springs*
- 13 Appliance Color Chart
- 14 Splints, Mouthguard and Pressure Formed
- 15 Perio Protect™ – *Your solution to healthy teeth and gums*
- 16 Laser Welding Certified
- 17 3D Imaging
- 18 Inman Aligner Clinical Manual
- 19 Education
- 20 Caring for Appliances and Cleaning Kits
- 21 Financial Requirements
- 22 Turn-Around Schedule & Shipping Forms
- 23 Shipping Form
- 24 Rx Forms

Inman Aligner Family of Appliances

The Inman Aligner Family of Appliances was developed in 2000 and earned the first patent in 2002. It is our flagship appliance and has given Inman Orthodontic Laboratories worldwide recognition. We are very proud of the success of the appliance line and we are sure there is a design that will fit your practice perfectly.



Comfortable!
Versatile!
Effective!



Inman Aligner Standard

The Inman Aligner is designed for maxillary and mandibular incisal correction. Lingual and labial components are used to correct anterior rotations and crowding. Used when force from the facial as well as the lingual is needed to push and pull (squeeze) the four incisors into alignment.



Clear Bow


The Clear Bow wire-free labial bow is a good option for patients concerned about aesthetics. Made from thermal-formed clear material the Clear Bow is almost invisible while being worn.

Inman Aligner Appliances



Inman Aligner Wraparound

This appliance is ideal when there is no room to cross the occlusion or when bridge work is involved. Fabricated with lingual and labial components to correct rotations and crowding.



Inman "DePaul" Habit Expansion

This design allows transverse expansion as well as anterior retraction as the habit is broken. It utilizes a unique trombone-like component to support the Bluegrass roller.



Inman Aligner Retraction Only

Used for retracting teeth to close spacing. The appliance may also be used to reposition teeth into a better position before cosmetic dentistry to reduce the amount of enamel reduction (preps) needed.



Inman Aligner Crossbite

Eliminates the need for finger springs, screws or adjustments. This comfortable appliance usually will correct anterior crossbites in less than three weeks.



Inman Aligner Habit

A great appliance used to close anterior spacing as tongue and finger habits are broken. When the teeth are moved into their ideal position the appliance becomes passive.



Inman Aligner "D"

This appliance consists of an exclusive pivoting anterior screw component that is connected to a palatal acrylic body and the Inman Aligner Retractable Bow. It is designed for anterior spacing cases in which incisor retraction and mesial bodily movement are required.



Invisible Inman Aligner


With this appliance, the lingual active component is fabricated to an invisible retainer. The Invisible Inman Aligner is used when lingual to labial movement is needed due to lingual displaced anterior teeth and aesthetics is a concern.



Inman Expansion Appliance

Lingual and labial coil components to correct anterior crowding and rotations. The expansion midline screw is added to gain space. This appliance will gain up to 3mm of arch expansion or upright the lower posteriors to create extra room to align the anterior teeth.

Inman Aligner – Adjusting the Appliance



Lingual Assembly Engaged



Lingual Assembly Removed

Lingual Component

The lingual component of the Inman Aligner is the real “magic” of the appliance. It creates force parallel to the occlusal plane allowing bodily movement of teeth, not just tipping. It creates a steady but gentle force by compressing NITI coil spring between the lingual acrylic pad and the stainless steel tubes.

The lingual component may be removed from the main body of the appliance by gently pulling the assembly toward the anterior.

Once the lingual assembly is fully disengaged, remove the original NITI coil spring. Replace the NITI coil spring with a heavier NITI coil spring to increase the force or a lighter NITI coil spring to decrease the force.

It is standard for us to assemble new Inman Aligners with **.010 X .030**.

A variety of wire gauges can be used to help achieve your result.

.010 X .030 Standard Forces
.012 X .030 Heavier Forces
.014 X .030 Heaviest Forces
.009 X .030 Reduce Forces

Labial Component

Because the force from the labial component is so direct it rarely needs to be adjusted, but can be easily adjusted if necessary.

If an adjustment bead was used in the fabrication process, grasp the helix of the terminal end of the wire passing through the buccal tube with bird beak pliers. Wind the helix in a clockwise direction to increase the force, counter clockwise to reduce the force.

On appliances without adjustment beads, cut the terminal end of the wire passing through the buccal tubes. Slide off the coil spring and replace it with a heavier NITI coil spring to increase the force or a lighter NITI coil spring to decrease the force.

A variety of wire gauges can be used to help achieve your result.

.010 X .030 Standard Forces
.012 X .030 Heavier Forces
.014 X .030 Heaviest Forces
.009 X .030 Reduce Forces



Labial component with stop.



Bend wire to compress coil.



Clasp Adjustment

Whenever possible the 2nd Bicuspid are clasped with an Adams Clasp with a laser welded tube. The position in the arch and the shape of the Bicuspid make them ideal for maximum retention. Sometimes due to a tight occlusion or crown and bridge work we will move the Adams Clasps to the molars or use a distal molar C-Clasp and cantilever the tube forward to the Bicuspid area.

Adams Clasps are adjusted by grasping the Mesial U Loop with Jarabak pliers. Twist the U Loop toward the cervical margin of the tooth to tighten the clasp. Repeat this to the Distal U loop to complete the adjustment process.

The Adams Clasps can become loose and deformed by the force transmitted to the clasps as the appliance is being removed. Over time the wires crossing the occlusion will need to be adjusted back down reversing the deformation, simple finger pressure is often all that is needed.

C- Clasps can also become distorted by the force transmitted to the clasps during removal. Simple finger pressure will return the clasp to the proper occlusal / gingival position. Using a pair of 3 prong pliers to bend the free end of the clasp tighter into the embrasure will increase the overall retention.



The Inman Aligner is sent for each patient with an Instruction and Care Pamphlet and an Appliance Case

The Requirements for Fabrication



- **Accurate maxillary and mandibular stone models. We need two models of each arch that will be treated with the Inman Aligner, or we can duplicate your models for an additional charge.**
- **A wax bite, registering proper articulation of teeth.**
- **A prescription form that indicates the design and teeth to be repositioned.**
- **Interproximal reduction: — mm.**
- **Adams clasp placement (standard 2nd bicuspids)**
- **Auxiliaries: coffin wire, tongue crib, tongue pearl.**

Inman Power Component

We designed the **Inman Power Component (IPC)** due to our frustration with the size and function of other components commonly used to compress open coil springs to generate force.

The concept and function of the **IPC** is simplistic beauty and allows adjustment by the patient or patient's guardian, making it a very efficient choice of creating force.



Multi-Functional!
Easy to Adjust!
Guaranteed to Never Slip!



Slotted Guide Bar Ni-Ti Open Coil Spring Sliding Lock


- The IPC is low profile and extremely comfortable compared to other components available.
- The lock is unidirectional and won't back up
- The IPC may be used for many types of movement, from simple space regainers to advanced 3-way arch development appliances.

Usage:

The infinite number of uses of this simple, dependable component will develop over time and we encourage you to use your imagination and knowledge to help us discover all of the **IPC's** potential. Inman Orthodontic Laboratory will work with you on any ideas you may have for the **IPC**. We will reward you for documented case results using your designs.

The Force:

The **IPC** compresses coil spring with an easy to adjust unidirectional sliding lock. The lock can be advanced 1mm "click" at a time for accurate, predictable activation.



Activating the IPC Appliances



The **IPC** is easily adjusted by advancing the sliding lock using our special **IPC** tool. A Ligature Director may also be used.

DO NOT ACTIVATE THIS APPLIANCE PRIOR TO INSERTION.

It is important to remember that the sliding lock is unidirectional and cannot be reversed. We use constant force Nickel Titanium Open Coil in all the appliances we fabricate; this allows you to have greater control of the tooth movement. We design the appliances with physical stops to control the amount of correction and eliminate the fear of over correction.


To reduce the force due to over-activation of the sliding lock, grasp one end of the open coil NiTi spring with Jarabak pliers and gently pull it. As the coil unwinds to the desired length, cut the coil flush with the component.

IPC Appliances



E-Arch

For controlled expansion of either arch, the IPC/Arnold Expander is a great choice. Insert it completely passive then adjust the sliding lock to generate the amount of force you desire. This is a proven design with the advantage of the control the IPC offers.



IPC Distalizer

We feel the IPC Distalizer is the easiest Distalizer to adjust. You can rely on 1mm of distalization per month and you can convert the appliance chairside from an active Distalizer to passive holding arch. To do so simply paint the IPC where it engages the tube with flowable composite and light cure, then cut and smooth the wire from the Nance button to the bicuspid band or rest.



IPC 3-Way

Our 3-way appliance uses a unique trombone mechanism that eliminates any cutting or modifications after expansion is achieved. Once you are ready to distalize simply slide the lock to compress the coil spring.

Arch Development

Developing arch form through the use of expansion is common and effective. Inman Orthodontic Laboratories can help you with our complete line of expansion appliances. We fit bands, add bite planes and add other auxiliaries to help you achieve great results with little time and hassle.

Please see the Inman and IPC Family of Appliances for more Arch Development designs.

Appliances



Removable 3-Way

For transverse expansion as well as A/P development. This appliance can be fabricated using one 3-way screw (for minor movements) or 3 individual screws for large corrections.



Hyrax

Rapid Expander for midline suture separation. It can be fabricated using many different size and style screws. It may be made using two bands and two rests (suggested) or four bands.



Schwartz Appliance

Removable option for expanding the upper and lower arch. The appliance can be fabricated with posterior bite blocks and many auxiliaries such as finger springs and habit controls.



Quad Helix

Highly adjustable expander that can be controlled by the clinician to the needs of the specific patient allowing differential expansion of the anterior and posterior.



Fan Screw Hyrax

Hyrax type screw to expand the anterior mainly and the posterior minimally.



Compact Hyrax

A compact and comfortable design utilizing 2 arms instead of 4.



Haas

This appliance adds palatal acrylic and is both tooth (bands) and tissue supported.



W-Arch

Similar to a Quad but without the helixes, may be used on either arch.



Bonded Hyrax

This appliance replaces the bands of a standard Hyrax and instead relies on bonding the occlusal acrylic to the posterior teeth. This design is often used in High Angle cases.



Flea

Fixed Lingual Expansion Appliance — a sturdy and reliable lower expansion design.



Fixed Sagittal

A great appliance to gain some A/P space as well as round out the arch. Adjustments are made by opening up the Omega Loops and generating anterior/posterior force.

*For more information,
call 1(800) 289.0118*

Whether you need to hold space, regain lost space or break habits, our metal appliances are carefully crafted and most connections are laser welded.

Holding Appliances



Space Regainer
An effective appliance for regaining space.



Band and Loop Space Maintainer
This appliance is used for maintaining space.



Lingual Arch
Used for the upper and lower arch, with or without the Omega loops. Used to maintain arch form.



Nance Holding Arch
This appliance keeps the molar position in developing dentition.



Trans Palatal Arch (TPA)
Keeps the molar position as well as holds expansion, the Omega loop may be activated to rotate the molars.

Habit Break Appliances



Tongue Crib
A great appliance for discouraging thumb sucking by blocking the thumb from the palate. This design may also be used for tongue thrust habits.



Blue Grass Roller
This is a counter conditioning appliance that eliminates the harmful thumb sucking habit by teaching the patient to play with the roller using their tongue. It is instant and very effective as the patients almost always find spinning the roller more satisfying than sucking their thumb.

These appliances are great for final retention. Most can be modified to aid in minor corrections. There is a variety of creative acrylic colors and patterns that patients will enjoy.


Retention Appliances



Spring Aligner
For minor tooth corrections of the centrals and laterals, this appliance may also be used as the final retainer.



Mini Spring Retainer
For very minor corrections of the centrals and lateral teeth as well as a final retainer. The addition of lingual clasps will reduce the chance of swallowing or aspiration, also they will help make the appliance harder for the patient to lose.



Hawley
Labial Bow cuspid to cuspid and an acrylic palate for retaining teeth after orthodontic treatment. Clasps, finger springs, labial acrylic and bite planes are some of the options available.



Retainers with Pontics
We can add a tooth or teeth to almost any retainer we fabricate. We need an accurate impression or model and a proper denture tooth (not ceramic) shade.

Wraparound Appliances



Wraparound
Labial bow from distal molar to distal molar eliminating any wires crossing the occlusion. This design allows settling of the occlusion and is a great choice for final retention. Clasps, finger springs, labial acrylic and bite planes are some of the options available.




Wraparound with Clasps
For better retention we can laser weld the labial bow to Adams clasps or C-clasp. Clasps, finger springs, labial acrylic and bite planes are some of the options available.



Wraparound with Support Wires
To keep the bow more stable we can add light wires between the laterals and cuspids or between the cuspids and bicuspid. This wire is tightly adapted and will not interfere with the occlusion.

3-3 Lingual Appliances



3-3 with Helix



3-3 with Composite Pads



3-3 with Mesh Pads

These appliances may be fabricated for the lower and the upper when occlusion allows. Since they are bonded, it eliminates compliance and all are good options for final retention.

We can add many different types of auxiliaries to any appliance. From extra clasps for increased retention to finger springs to move a tooth or teeth. We can also incorporate labial acrylic to the bow to eliminate any chance of rotation or screws for expansion. Here are examples of some of the auxiliary options we offer.

Retention



Adams Clasp



Arrow Clasp



Ball Clasp



Adams Clasp laser welded to wire



C-Clasp laser welded to wire



C-Clasp laser welded to labial bow

Bite Blocks



Bite Blocks and Screws



Screws

Finger Springs



Distal Kicker Spring




Closing Spring



Z Spring



Sweep Springs



Helix Spring

We make wearing appliances fun and exciting for your patients. Choose from a variety of wild colors and custom designs.



Our mouthguards, splints and pressured formed appliances are strong, reliable and offer unparalleled patient comfort.

Splints



Flat Plane Splint

Used for TMJ problems as well as to protect dental work or natural dentition from bruxing.



Flat Plane with Anterior Guidance

This splint is designed so that as the patient moves through excursions the posterior teeth lose their contact (disocclude) with the splint.



Deprogrammer

This appliance only allows contact with the lower centrals and laterals not allowing the patient to clench and cause muscle disharmony. This is a very effective and simple appliance and is very well accepted by patients.

Pressure Formed Appliances



Athletic Mouthguard

All your patients involved in sports should be advised to wear a custom made, fitted mouthguard to prevent tooth, neck and jaw damage as well as concussions and other serious injuries. We can fabricate mouthguards in different thickness and a variety of colors.



Bleaching Trays

We use BioStar Positive Pressure Forming machines to fabricate well adapted and tightly sealed bleaching trays.



Invisible Retainers / Essex Retainers

Our standard is .040 or 1mm before being formed but we can use any thickness you prefer. Teeth can be added for support of spaces.


Give Your Patients an Alternative in Treating Gum Disease.

What is Perio Protect™

Perio Protect™—Your solution to healthy teeth and gums. Perio Protect™ is a revolutionary new method to treat periodontal disease using an FDA-cleared medical device to deliver dentist-prescribed medications directly to the source of the infection. These medications address the pathogens that cause gum disease.

The medical device is a set of custom trays prescribed by the doctor and fabricated in our laboratory. The trays offer varying thickness and differing extensions to support a special patented seal that directs antimicrobial agents of the dentist's choosing into the gingival sulcus, and the method maintains the medications for a sufficient frequency and duration of time to allow the medications to control the pathogens causing the infection.

As healing occurs, patients alter the frequency and duration that they wear their trays.



For more information on this great product call Inman Orthodontic Laboratories 1(800) 289.0118

For more information, call 1(800) 289.0118



Before treatment with Perio Protect™ Trays



After treatment with Perio Protect™ Trays

Laser Welding Certified

Because we always strive to fabricate the best possible appliances, Inman Orthodontic Laboratories laser welds all of our metal to create the strongest appliance possible, without compromising the integrity of the wire.

Inman Orthodontic Laboratories Utilizes Laser Welding



The Difference...



Laser Weld

Laser Welding is a process in which light energy is used to join parent metals to themselves.

- This process fuses the parent metal on a molecular level resulting in a finished product that is one alloy.
- The heat used in this process is so localized that it results in a seamless, undetectable work zone that is not discolored or compromised in any way.
- Eliminates porosity.

Solder Joint

Soldering is a capillary fill system where solder is heated with a flame until it flows, joining the two metals together.

- Solder is designed to melt at a lower temperature than the parent metal, causing weakness.
- Solder is a different and much weaker alloy than metal.
- The heat used for soldering is very high often resulting in unwanted tempering of the wires near and in the soldered area.

3D Imaging


Inman Orthodontic Laboratory Digital Study Models with advanced diagnostic tools.



3-axis, two camera laser 3D scanner



Original Model



3D Scan

3D Imaging

- Highly accurate
- Full coverage — scans all undercuts automatically using 3-axis motion system
- Easily scans the antagonist

Advantages of 3D Image

- Color Distance Mapping
- Occlusion analysis
- Measuring and analysis tools
- Tooth width analysis
- Bolton analysis
- Space analysis
- Ideal arch analysis/available space

NEW! 3-D Digital Printing

Inman Orthodontic Laboratory continues to invest in new technology to better serve our Doctors.

We can create 3D models from digital dental impressions or digital scans of dental models.

Our *Objet Rapid Prototype* machine can "print" accurate 3 dimensional dental models that can be used to fabricate appliances.



We can also do diagnostic set-ups and print the models in various stages of correction allowing treatment planning options.



From a digital file to a 3 dimensional dental model


To learn more about any of these services, please call 1(800) 289.0118

The 160 page Clinical Manual for the Inman Aligner™

This hard cover manual explains how to incorporate the Inman Aligner into your dental practice by illustrating proper Case Selection techniques and Arch Analysis.

With the use of this manual, you will have the confidence to add a new service to offer your patients and increase your bottom line.

Everything is covered from Anchors to Spacewize™.



The manual includes:

- The Workings of the Inman Aligner
- Model Evaluation and Arch Analysis
- How to do Interproximal Reduction
- Correcting Crowding
- Case Studies
- Clinical Cases
- Diagnostic and Treatment Plan Review
- and much more!



Before treatment with Inman Aligner™



After treatment with Inman Aligner™

The Inman Aligner™ Hands-On Seminars

These thoroughly comprehensive Seminars are offered periodically at various locations throughout the year.

Learning Objectives:

- Diagnosis and case selection
- Movement potential
- Arch evaluation
- Fitting and adjusting an aligner
- Interproximal reduction
- Composite anchor placement
- Retention technique
- Restorative pre-alignment
- Ethical considerations

"Incredibly comprehensive, simply brilliant. Truly inspirational!"
Dr. Sverker Toreskog, Internally renowned speaker



These courses are reserved to 40 dentists and 20 staff.

The Inman Aligner™ On-Line Course



If you can't attend a seminar, take the Inman Aligner on-line course.

www.straight-talks.com

The Inman Aligner on-line course includes:

- Concept and Philosophy
- Key Learning Concepts
- Case Studies
- Animation and Videos
- Techniques
- Learning-Aid Quiz
- ... and much more

The Inman Aligner™ On-Line Support

With the certification, you will have access to our on-line support and tools.

- Comprehensive articles and information
- Post your cases and discuss with other dentists
- Utilize our Spacewize™ calculator- instantly accesses cases at your consultation



Articles and Videos



Spacewize™ Calculator



Visit www.straight-talks.com and begin your Inman Aligner course study program today.

Keep Your Appliances Clean and Fresh!

Retainer Brite cleans your Hawley Retainers, Inman Aligners, Mouthguards, Splints, and most other removable Dental Appliances.

- Helps remove plaque and tartar
- Kills germs and bacteria
- Cleans and brightens
- Keeps appliances fresh



Call to order
1(800) 289.0118

Retainer Brite Cleaning Tablets and Sonic Cleaner

includes the following:

- Sonic Cleaner
- 1 year supply of tablets



Retainer Brite Cleaning Tablets

- 36 count – Case of 24
- 92 count – Case of 12

Your account is opened with your first order. It's that simple!

You will receive a welcome packet including Rx forms, mailing labels, and information regarding ordering, scheduling, and procedures. As a courtesy, we will send you a monthly statement summarizing your activity.

Payment upon receipt of your invoice is greatly appreciated.

Your account is due upon receipt of your statement and credit card payments are accepted through the 15th of the month. We suggest sending the right hand portion of your statement along with your payment or writing the invoice number on your check.



Experience the difference...
your patients will appreciate the results!

PAYMENT METHODS

AUTO-BILLING: your credit card will be automatically debited from your financial institution and your statement paid on a date set with our office manager.

CREDIT CARD PAYMENTS: We accept Visa, MasterCard, Discover, and American Express. To pay by credit card, please write the account number, expiration date, three digit security code located on the back of your MasterCard, Visa, or Discover. Your American Express card has a 4 digit code located above the account number.

Invoices over 30 days will incur a 1.5% late charge (19.56% per annum), which will be posted to your next statement. Any account over 90 days will be processed through a collection agency and your account will be placed on C.O.D. status. These procedures are necessary as delinquent accounts place a severe handicap on the office overhead and cash flow, which eventually results in higher costs and increased fees to our clients.

*Please call to speak to a friendly member of our staff!
(954) 340-8477 or
1(800) 289-0118.*

We are always here to assist you!

Convenient, Reliable Shipping

Local Accounts

- Count 5 working days for appliances, beginning with the day after the package leaves your office.
- Local courier service is available. To schedule a pick-up for the same day you must call before 1 p.m. To schedule a pick-up for the next day you must call before 5 p.m.

Accounts Outside of our South Florida Area

- Count 9 working days for appliances beginning with the day after the package leaves your office.
- UPS 2nd day delivery is available to all out of area clients.

Please use the production turn-around chart and the work calendar to assist with your scheduling.

APPLIANCES	PRODUCTION TIME
Retainer	3 business days
Metal Appliances	3 business days
Splints	4 business days
Herbst Appliance	9 business days
Functionals	5 business days
Inman Aligner	7 business days
Perio Protect Trays	5 business days




Please call 1(800) 289.0118 to order your mailing packet of Rx forms and mailing labels.



Mailing Label

*Please copy and use the appropriate Rx form and/or mailing label for your first case.
Please call 1(800) 289.0118 to replenish your mailing packet of Rx forms and mailing labels.*



Rx Forms

Rx forms can be copied, filled out and sent with your impressions or models.
Call us at 1(800) 289.0118, if you need any mailing supplies.



Experience the difference...
your patients will appreciate
the results!



DL #0001147

Dr. Name _____ Patient Name _____
Address _____ Date Sent _____
City _____ Date Due _____
State/Zip _____ License # _____

INMAN ALIGNER RX

Reset

21	12
21	12

Only on very rare occasions
can we change the position
of the cuspsids.

- Strip as Necessary
- Minimal Stripping
- Inman Cross Bite
- Inman Habit Fixed
- Inman Standard
- Inman Fixed Bow
- Wraparound Design
- Retraction Only
- Expansion Screw
- Adjustment Beads (Pair)


SUPPLIES

- Spool of NITI Coil
- 010 x 030 Patient Information Pack
- 012 x 030 Mailing Supplies
- Appliance Cleaning Kit

COLORS:

Upper _____ Lower _____

Special Notes: _____



Experience the difference...
Your patients will appreciate
the results!




DL #0001147

Dr. Name _____ Patient Name _____
Address _____ Date Sent _____
City _____ Date Due _____
State/Zip _____ License # _____

RX

- Hawley
- Wraparound
- Modified Spring Ret. reset
- Schwartz
- Sagittal
- 3-Way
- Splint
- Quad
- Lingual Arch
- Habit Appliance
- Hyrax
- Space Maintainer
- Pedo Partial Shade
- Fixed Sagittal
- W-Arch
- Bionator-Twin Block-OC1-OC2
- Herbst-Crowns
- Bands
- Tubes
- Ball Clasps (Between)
- Adams Clasps (Teeth)
- C-Clasps (Teeth)
- Elbow (Between)
- Soldered C-Clasps (Teeth)
- Springs (Teeth)
- Springs (Teeth)



COLORS:

- Regular Upper _____ Lower _____
- Tropical Upper _____ Lower _____
- Neon Upper _____ Lower _____
- Electric Upper _____ Lower _____
- Graphic Upper _____ Lower _____
- Custom Upper _____ Lower _____

Special Notes: _____

LAB USE	
Tech. Int.	
Wire	_____
Acyl	_____
Trim	_____
Pol	_____
Pak	_____

DENTIST SIGNATURE _____

3953 NW 126th Avenue • Coral Springs, FL 33065 • (954) 340-8477
A Certified Dental Laboratory

CALL TOLL FREE 1-800-289-0118

We need two models of the arch involved and a counter model or pvs impressions. If they are not provided a duplication fee will be added.

DENTIST SIGNATURE _____

3953 NW 126th Avenue • Coral Springs, FL 33065 • (954) 340-8477
A Certified Dental Laboratory

CALL TOLL FREE 1-800-289-0118